

SAPPHIRE INFUSION PLATFORM US CATALOGUE 2019



SAPPHIRE
by Eitan Group

CONTENTS

■	<u>INFUSION KITS</u>	4
■	<u>ACCESSORIES</u>	6
■	<u>ADMINISTRATION SETS</u>	10
■	<u>PC APPLICATIONS</u>	13

Sapphire™ Multi-therapy infusion pump kit



Includes:

- Sapphire™ Multi-therapy pump
- Power supply
- Mini cradle
- Mini cradle knob key
- User manual

P/N 17000-028-0072

Sapphire™ Epidural infusion pump kit



Includes:

- Sapphire™ Epidural pump
- Power supply
- Bolus handle
- Mini cradle with splitter
- Mini cradle knob key
- User manual

P/N 17000-031-0035

Sapphire™
Mini Cradle



Sapphire™
Mini Cradle

With Splitter

P/N 17000-026-0007

Without Splitter

P/N 15042-000-0010

Mini Cradle with IPS

P/N 15042-043-0025

Sapphire™
Multi Pump Mounting System



P/N 17000-020-0008

Sapphire™
Integrated Power Supply



P/N 15072-000-0010

Sapphire™
Power Supply



P/N 15072-000-0007

Sapphire™
External Battery Pack (EBP)



P/N 15079-111-0001

Sapphire™
AC Power Cord - Type B



P/N 15070-112-0005

Sapphire™
Power Cord - compatible
with IPS (US and Canada)



P/N 15171-112-0002

ACCESSORIES | BAGS

Sapphire™
Infusion Pouch 500 mL



P/N 05080-605-0001

Sapphire™
Travel Case



P/N 15020-010-0002

Sapphire™ Homecare Large Backpack



P/N 15088-000-0001

ACCESSORIES | PAIN MANAGEMENT

Sapphire™
PCA lockbox 250mL



P/N 15063-000-0001

Sapphire™
PCA lockbox 100mL






P/N 15061-131-0004

Sapphire™ Bolus handle






P/N 15071-000-0005




ADMINISTRATION SETS | STANDARD SETS

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/ BOX
AP 403-01	12003-000-0013	Microbore infusion set with non-vented spike	290	~6 mL	30
					
AP 409-01	12005-000-0003	Primary infusion set with vented/non-vented drip chamber, back check valve and 2 needleless y-site	277	~23 mL	30
					
AP 423-01	12003-000-0023	Microbore with non-vented spike and a back check valve	292	~7 mL	30
					
AP 433-01	12000-000-0044	Macrobore Infusion set with vented/non-vented drip chamber and needleless y-site	286	~22 mL	20
					

ADMINISTRATION SETS | FILTER SETS

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/ BOX
AP 204-01	12003-000-0070	Microbore with non-vented spike, 0.2µm filter and back check valve	295	~8 mL	30
					
AP 206-01	12000-000-0022	Macrobore with non vented spike, 1.2µm filter and needleless Y-site	295	~23 mL	20
					
AP 214-01	12004-000-0002	Epidural (yellow marked) microbore with vented/non-vented spike and 0.2µm filter with non-vented spike, side arm with female luer lock and back check valve	295	~ 8 mL	30
					
AP 240-01	12003-000-0069	Microbore with non-vented spike, 1.2µm filter and back check valve	295	~ 8 mL	20
					

ADMINISTRATION SETS | SPECIALTY SETS

NAME	P/N	DESCRIPTION	LENGTH (CM)	PRIME VOLUME	UNITS/ BOX
AP 133-01	12006-000-0010	Blood set with 2 non-vented spikes,200um filter and needleless y-site	292	30 mL	20
					
AP 402-01	12004-000-0004	Epidural (yellow marked) microbore with non-vented spike	290	~6 mL	30
					
AP 410-01	12005-000-0007	Secondary infusion set with vented-non vented drip chamber	71	~9 mL	30
					
AP 424-01	12003-000-0024	Microbore with non-vented spike, side arm with female luer lock and back check valve	294	~7 mL	30
					
AP 504-01	12000-000-0045	UV-Resistant macrobore with vented/non-vented drip chamber and needleless Y-site	287	~26 mL	20
					

PC APPLICATIONS

The Sapphire modern, software based platform expands beyond the pump to the PC applications suite that complements the pumps;



PUMP LOADER

Enables upgrading or reloading software to the Sapphire pump via a wire connection.



DRUG LIBRARY

Enables to set configuration settings, drug specific limit and more. When complete, the library can be uploaded to multiple Sapphire pumps.



EVENT VIEWER

Enables downloading and viewing the event log of the Sapphire infusion pump on a computer.



ANNUAL CERTIFICATION KIT



(SW + disposable kit [15075-020-0001])
Enables to run performance test (air detection, pressure detection and accuracy). Results are compared to manufacturing database per pump.



PC APPLICATIONS | ACCESSORIES

PC APPLICATION	ITEM	P/N	PICTURE
Annual Certification Test	Annual certification kit	15075-020-0001	
All	Communication cable	05020-110-0213	
All	RS232 to USB 4 port adapter	15077-000-0001	

ACCESSORIES | SPARE PARTS

SPARE PART FOR	ITEM	P/N	PICTURE
Sapphire pump	Li-Ion battery	05020-160-0001	
Mini cradle	Splitter	15043-000-0005	
Mini cradle	connection knob flange Key	05040-380-0019	
Mini cradle	connection knob key	05040-380-0018	
PCA lockbox 100	Key	05070-190-0015	
PCA lockbox 250	Key	10010-122-0046	
PCA lockbox 250	Release clip (eliminates lockbox locking onto cradle)	15063-121-0029	
PCA Lockbox 500	Key	15062-132-0001	

REF 10050-003-0009 | Rev.00 / 06.2019

Q Core Medical Ltd.
29 Yad Haruzim St., P.O. Box 8639
Netanya 4250529, Israel
Tel: +972.73.2388888
Fax: +972.73.2388800
Email: support@Qcore.com
www.qcore.com



SAPPHIRE
by Eitan Group